


## X-Ray Fluorescence Analysis of Chlorine in Cement and Hot Meal Using Borate Fusion and the VITRIOX® electric

### Introduction

The analysis of materials containing chlorine with XRF is a continuous challenge, because chlorine displays high volatility during borate fusion. This leads to poor repeatability with conventional fusion systems.



**Fig. 1: Electrical Fusion Machine with 4 Stations for XRF and ICP.**

## X-Ray Fluorescence Analysis of Chlorine in Cement and Hot Meal Using Borate Fusion and the VITRIOX® electric

The method presented here takes advantage of the capabilities of the new electrical fusion machine VITRIOX electric from FLUXANA [1]:

- Precise temperature control
- High precision
- Fusion with covers



**Fig.2: Cup for electrical fusion machine with removable cover.**

### Procedure

Sample preparation of the sample, which had been dried at 105 °C, was conducted using borate fusion. The ratio of sample to flux was 1:8. In this way, it was possible to achieve high sensitivity for chlorine.

### Sample Preparation


Cement, hot meal dried            1 g

Flux FX-X65\*                            8 g

\*66% lithium tetraborate + 34% lithium metaborate

### X-Ray Fluorescence Analysis of Chlorine in Cement and Hot Meal Using Borate Fusion and the VITRIOX® electric

#### 1. Example: Calibration of Chlorine in hot meal



**Fig. 3: Calibration for chlorine in hot meal, calibration error RMS = 0.02%. Reference values were performed with titration.**

The calibration samples were real samples provided by a customer analyzed by wet chemistry. They were fused using the electrical fusion machine from FLUXANA. Duplicates of all samples were produced. The calibration error achieved for chlorine was 0.02%.

## X-Ray Fluorescence Analysis of Chlorine in Cement and Hot Meal Using Borate Fusion and the VITRIOX® electric

### 2. Example: Determination of chlorine in cement

Table 1 shows the results for repeat preparations with the VITRIOX® Electric of cement samples containing low chlorine concentrations.

**Table 1: Determination of chlorine in cement using fused bead method.**

Preparation#	Chlorine mass%
1	0,037
2	0,036
3	0,038
4	0,035
5	0,036
6	0,035
7	0,037
8	0,036
9	0,034
10	0,035
11	0,034
12	0,041
<b>Mean</b>	0,036
<b>Std.dev.</b>	<b>0,002</b>

### Summary

The results presented here clearly confirm that borate beads can be produced with highest precision using FLUXANA's new electrical fusion machine VITRIOX® electric. Volatile elements, such as chlorine, for which the precision depends on the temperature stability of the fusion machine, can be satisfactorily analyzed.

### References

[1] Rainer Schramm, X-Ray Fluorescence Analysis: Practical and Easy - 2nd edition, FLUXANA (2017).